

PRV

PATENT- OCH REGISTRERINGSVERKET

Patentavdelningen

**Intyg
Certificate**

Härmed intygas att bifogade kopior överensstämmer med de handlingar som ursprungligen ingivits till Patent- och registreringsverket i nedannämnda ansökan.

Ansökan ingavs ursprungligen på engelska.

This is to certify that the annexed is a true copy of the documents as originally filed with the Patent- and Registration Office in connection with the following patent application.

The application was originally filed in English.


(71) Sökande AstraZeneca AB, Södertälje SE
Applicant (s)

(21) Patentansökningsnummer 0303030-1
Patent application number

(86) Ingivningsdatum 2003-11-14
Date of filing

Stockholm, 2004-11-15

För Patent- och registreringsverket
For the Patent- and Registration Office


Gunilla Larsson

Avgift
Fee

**PATENT- OCH
REGISTRERINGSVERKET
SWEDEN**

Postadress/Adress
Box 5055
S-102 42 STOCKHOLM

Telefon/Phone
+46 8 782 25 00
Vx 08-782 25 00

Telex
17978
PATOREG S

Telefax
+46 8 666 02 86
08-666 02 86

Summary of Invention

The invention relates to the removal of the contents of an MDI (metered dose inhaler) canister without causing the turbulence associated with this activity when conducted by the normal means. In doing so, the device allows the accurate determination of the amount of drug adhering to the inner side of the canister and/or the valve parts separately and in addition to the drug content of the expelled portion of the canister.

Background to the Invention

It is well known in the field of metered dose inhalers that some drugs adhere to the sides of the canister and onto valve parts in contact with the formulation. The formulation is usually in the form of a suspension of drug and/or excipients in a volatile propellant. Information on the amount of drug adhering to internal surfaces of the manufacturing equipment and in the pMDI is essential in the development of the product.

The most favoured way of alleviating this problem is to coat the inside of the can with a low energy coating such as those described in WO01/51222 and EP642992.

Drug adhesion is a significant problem in that it reduces the availability of the drug to the patient. Measurement of the parameter of potential drug adhesion allows different coatings to be explored and makes calculations of the drug overage to be accurately ascertained. The manufacturing overage is calculated using the concentration of the drug in the manufacturing vessel. An overage is the extra drug which has to be added to the formulation during manufacture to account for that which may adhere to the manufacturing vessel and the inside of the canister over the time it is stored prior to reaching the patient. Regulatory bodies generally require an overage of less than 10%. Accurate determination of loss of drug during manufacture and on the wall

of the canister are therefore essential.

Also, the tendency for the pMDI to be stored inverted, i.e. valve down, inside the actuator exacerbates the problem as the formulation is in contact with the valve for a significant time.

The potential risks of a high can overage comes from the possibility of reversibility, where the excess drug falls back into suspension and results in a higher than required dose of drug delivered to the patient. The dose uniformity testing carried out during batch testing may determine this but would result in batch failure.

It is imperative, therefore, that the amount of drug adhesion (if any) is determined at an early stage of the MDI development. The normal procedure of cooling the can, piercing it and pouring out the contents is not accurate as the cooling itself can cause deposition (as a result of propellant evaporation) or cracking of the deposit, (in the latter, causing it to drop into the main body of the formulation), each in turn giving a falsely high or low result.

Piercing the can and allowing the contents to escape under pressure, is not a viable alternative either as the drug may be deposited as the propellant evaporates or is disturbed and removed during the uncontrolled exit of the propellant. Therefore, it can be seen that current methods for determining the amount of can deposition are not reproducible and are not suitable for low dose products as small losses in the method become unacceptable when small amounts are to be recovered i.e. micrograms.

There are several can piercing applications in the prior art, for example, US patent number 3828,976. However the majority are for emptying of the can per se, either for capture of the propellant or recycling of the can material and are not intended to be a vehicle by which the can and contents are analysed.

US patent US 6393,900 details an apparatus for emptying the contents of an MDI and capturing the contents for analysis. This is essentially an automated can contents

analysis and gives total can concentrations. The invention does not use a pressurised system and is not designed to measure drug adhesion inside the canister as it only takes into account the total of the can contents. The contents are forced out of the canister under its own pressure; there is no effort to control it only to contain it.

The invention detailed in this application allows a differential analysis of the can contents, including the proportion of drug left inside the canister and, (if the valve is detached before the empty canister is washed), a measurement of drug on can and valve separately.

Another patent, DE20203999, details an apparatus which pierces the canister and confines the contents for analysis but the canister is not pre-pressurised and the primary objective is to empty the can into a receptacle without having to cool it first and to use an ion-selective electrode to determine a required parameter under pressure while the contents are held within the receiving vessel. There is no attempt to measure contents left inside the canister.

Therefore, it is the purpose of this invention to control the exit of the contents of the canister from said canister in a controlled manner such that the layer of drug adhering to the inside of the can and on the valve mechanism may be left behind unperturbed and therefore accurately measured.

Description of the Invention.

One form of the present invention may broadly be said to consist in a method of analysing the contents of a pressurised container comprising the steps of: enclosing said container in a pressure vessel; pressurising said pressure vessel with a non-reactive fluid; piercing said pressurised container within said pressure vessel; and analysing the content of said pressurised container when drawn off through said piecing.

Preferably the pressurised container is a canister container medicament.

Preferably said canister is a metered dose inhaler canister.

Preferably said non-reactive fluid is nitrogen.

The apparatus is a rig for piercing an MDI canister in order to empty it of its contents without disturbing any material adhered to the surface of the can.

If rapid emptying took place, i.e. as pressure is suddenly released, the can deposition would be disturbed. In this instance the can is pre-pressurised by being pierced in the side, initially by a cannula under pressure by nitrogen.

As the base is then pierced, the pressure inside the can is kept constant allowing the can contents to flow out of the base of the canister with control, thereby not disturbing any material deposited on the sides or on the valve.

Obstacles to be overcome in the design of the invention were:

Designing a vessel which could hold a pressurised canister securely enough to pierce it using external controls without a leak or loss of can contents into the holding vessel.

Minimising losses in the transfer line linking the piercer with the collection vessel.

To validate the wash efficiency to ensure full recovery of deposited and evacuated drug.

Devise a means of evacuating the collection vessel of propellant without losing drug content.

To pierce the canister and provide a pressurised supply to the inside of the can, prior to piercing at the base

To be able to determine accurately the drug left on the can wall separately from the rest of the can contents.

The following apparatus was devised to overcome these obstacles:

The device comprises several components:

- a) a pressuriser comprising a nitrogen supply, control gauge and connecting tube
- b) a collection vessel comprising three sections; (i) the main body and lid (ii) a three-way valve on the top (iii) a nitrile seal between the base and lid to contain and maintain pressure until venting is needed.
- c) A two compartment device which constitutes the main part of the invention comprising two halves joined together by three screws, with an internal void of suitable size to accommodate the pMDI canister. The pMDI fits in tightly and is sealed in place by three rubber O-rings. The 3 screws seal the apparatus when tightened.

The top section has a rotating bar attached at the side. When rotated the piercer enters the space occupied by the pMDI canister. The piercing action is similar to a needle. Once the can is pierced the nitrogen pressure is applied via the connecting tube. The bottom section has a piercer at the base and with the aid of a second rotating bar; this pierces the base of the canister and is then retracted.

The base has a hole at the bottom corner, which is so designed to channel the product via the transfer line under pressure from the nitrogen, into the collection vessel.

Method of Sample Collection to determine Drug adhesion.

To summarise, the method of collection is as follows:

Weigh the can and vessel

Seal can into piercer

Pierce side of can allowing nitrogen to enter at required pressure

Pierce base of can and retract
Open 3 way valves allowing product to flow into collection vessel
Leave nitrogen flowing for several minutes
Seal vessel and can piercer using 3 way valves
Reweigh can and vessel
Vent into dose unit (only on validation)
Wash vessel with solvent to dissolve drug and excipients
Collect and analyse drug remaining in the can and on valve.

The wash method will be determined by the product and would be determined as a normal part of validation of the recovery of the drug from the system and is well within the knowledge of a competent person skilled in analytical chemistry. The drug content assay would also constitute normal practice as part of the usual product development process and within the scope of a skilled person.

EXPERIMENTAL DETAILS

The formulation under test was a low concentration suspension formulation of formoterol fumarate dihydrate (FFD). The drug was suspended in a blend of propellants HFA 134A and HFA 227. The description of the formulation is contained in patent application WO03/63843.

The main concerns of the experimental validation were a) leakage b) recovery efficiency. Parameters affecting recovery were: efficiency of transfer along the transfer line linking the can piercing unit with the recovery vessel, the vessel washing procedure, potential losses on venting the collection chamber when venting into a dose collection vessel and subsequent washing of that vessel, and recovery of the deposits of drug inside the vented canister.

Placebo cans were used to check for carryover of FFD from one can to the next. Negligible amounts of FFD were found using the placebo cans showing that carryover of drug was not a concern and that the wash method was efficient.

Early trials resulted in low recovery due to leakage of the O-rings, which sealed the vessels with subsequent loss of drug. Care must be taken to align the O-rings correctly and seal the vessel. Any damaged O-rings must be replaced immediately.

To this end the pressure was lowered from 7 bar to 4 bar and the nitrogen left flowing for 1 minute into the piercing device prior to piercing the can. The effect was to reduce the FFD deposited in the piercer.

Summary of the Experimental Validation of the Device

1. Placebo cans showed that there was negligible drug deposited in the transfer line.
2. When the pressurised vessel was vented into a dose delivery unit on venting to test for loss of drug during the venting process, the resulting wash solution showed no traces of drug.
3. The wash procedure for the collection vessel used three washes. Negligible amounts of drug were found in the second and third washes, meaning a single wash could be used. (This, of course, may vary from drug to drug)
4. The can piercer was analysed for FFD deposits. Initially it was found to contain a significant amount of drug. However allowing the piercer device to pressurise for 1 minute prior to piercing and lowering the pressure (thereby reducing the pressure differential between can and device) eliminated this problem. To be certain of full recovery of the FFD, this was washed as well.

The validation and optimisation of the method resulted in an acceptable recovery of FFD, such that quantisation of can deposition could be determined reliably.

Results of cans analysed by the resulting method can be seen in Table 1 and the method is outlined below:

Weigh the can and vessel

Seal can into piercer as described previously

Leave nitrogen flowing at required pressure for 1 minute

Pierce side of can allowing nitrogen to enter at required pressure

Pierce base of can and retract

Open 3 way valves allowing product to flow into collection vessel

Leave nitrogen flowing for several minutes

Seal vessel and can piercer using 3 way valves

Reweigh can and vessel

Vent collection vessel

Wash vessel and can piercer with solvent to dissolve drug and excipients

Collect and analyse drug remaining in the can and on valve.

Calculations used in Tables 1 and 2

A = Residual drug in Can (% w/w)

B = Vessel Washings (% w/w)

C = Can Piercer Washings (% w/w)

D = A+B+C = Total FFD Recovered from can and piercer components (% w/w)

E = Total recovery = $D/F \times 100$ (%)

Where F is the expected total can content (from QA batch testing data). For uncoated cans using 134A, expected total can content was 0.0216 %w/w and for blend formulation in coated cans it was 0.0167 % w/w.

TABLE 1: results of cans tested using the system described above.

DRUG IN UNCOATED CAN IN HFA134A

Ca n	Residual drug in Can % w/w	Vessel Washings % w/w	Can Piercer Washings % w/w	Total FFD Recovered % w/w	Recovery %	Conditions
1	0.0086	0.009	0.0018	0.0194	90 %	4 bar pressure
2	0.0080	0.0101	0.0019	0.0200	93 %	4 bar pressure
3	0.0117	0.0075	0.0019	0.0211	98 %	7 bar pressure
4	0.0114	0.0086	0.0019	0.0219	101 %	7 bar pressure
5	0.0096	0.0097	0.0019	0.0212	98 %	4 bar pressure
6	0.0106	0.0039	0.004	0.0185	86 %	4 bar pressure
7	0.0092	0.0102	0.0023	0.0217	100 %	4 bar pressure
8	0.0098	0.0069	0.0027	0.0194	90 %	4 bar pressure
9	0.0096	0.0097	0.0015	0.0208	96 %	4 bar pressure

TABLE 2: Results

Drug in coated can in HFA 227/134A blend

Ca n	Residual drug in Can % w/w	Vessel Washings % w/w	Can Piercer Washings % w/w	Total FFD Recovered % w/w	Recovery %	Conditions
1	0.0037	0.0095	0.0028	0.016	96 %	4 bar pressure
2	0.0027	0.0094	0.0023	0.0144	86 %	4 bar pressure
3	0.0028	0.0056	0.0044	0.0128	77 %	4 bar pressure
4	0.0025	0.0128	0.0018	0.0171	102 %	4 bar pressure
5	0.0024	0.0071	0.0042	0.0137	82 %	4 bar pressure

The most noticeable difference in the two batches is between the uncoated and coated cans: the deposition in the uncoated can having on average 0.0100 w/w % deposition and the average on the coated can was 0.0028 w/w%, this represents a significant effect when using a coated can. I.e. the can deposition in the 134A/uncoated can is, on average, 45% compared to 17% in the blend/coated can, when compared with the total Formoterol content in the formulation.

It is entirely possible that the device and method would therefore distinguish between different formulation deposits using the same coating or the same formulation exposed to different coatings. The invention is not limited to the formulations or canisters described here and could easily be applied to other cans of differing size by altering the dimensions of the can holding vessel.

The conclusions which can be drawn from the results in Tables 1 and 2 is that:

The device is suitable for emptying can contents controllably such that any drug deposition on the canister wall is left undisturbed and can be analysed separately.

The method is accurate and reproducible and can be used to ascertain the differences between deposits of drug in different environments such as uncoated and coated surfaces or between formulations.

Separation of can and valve allows evaluation of deposition in each component.

What we claim is:

- 5 1. A method of analysing the contents of a pressurised container comprising the steps of: enclosing said container in a pressure vessel;
pressuring said pressure vessel with a non-reactive fluid;
piercing said pressurised container within said pressure vessel;
and analysing the content of said pressurised container when drawn off through
10 said piercing.
2. A method of analysing the contents of a pressurised container as claimed in claim 1
wherein the pressurised container is a canister container medicament.
- 15 3. A method of analysing the contents of a pressurised container as claimed in claim 1
or 2 wherein said canister is a metered dose inhaler canister.
- 20 4. A method of analysing the contents of a pressurised container as claimed in any one
of claims 1 to 3 wherein said non-reactive fluid is nitrogen.

Abstract

10

Document made available under the Patent Cooperation Treaty (PCT)

International application number: PCT/SE04/001642

International filing date: 11 November 2004 (11.11.2004)

Document type: Certified copy of priority document

Document details: Country/Office: SE
Number: 0303030-1
Filing date: 14 November 2003 (14.11.2003)

Date of receipt at the International Bureau: 26 November 2004 (26.11.2004)

Remark: Priority document submitted or transmitted to the International Bureau in compliance with Rule 17.1(a) or (b)



World Intellectual Property Organization (WIPO) - Geneva, Switzerland
Organisation Mondiale de la Propriété Intellectuelle (OMPI) - Genève, Suisse